


INTRODUCTION



- Particulate contamination from single-use systems (SUS) for pharma/biopharma production can potentially impose negative risks for yield, purity, and patient health¹.
- Effective particulate contamination monitoring and control are important to minimize the risks.
- Pharmaceutical products need to meet USP <788> and USP <790> criteria to get FDA approval for patient use.
- Entegris' Liquid Particle Counting (LPC) systems can measure particle count with size down to 0.5 μm .
- Entegris' Aramus fluoropolymer SUS is much cleaner than other SUS with different materials.
- Entegris is positioned uniquely to provide both SUS with minimal particulate contamination and instruments to monitor particulate level in SUS.

ARAMUS SUS PARTICULATE TESTING & BENCHMARKING



- Aramus SUS passes USP <788> - Particulate matter in injections
- PE bag of Brand B and Brand C higher counts with large standard deviation, some fails USP <788> (inset)

WHY PARTICULATES MATTER

Process performance – **yield**

Product quality and efficacy – **purity**


Patient safety – **health**

Harm = Impact of particulate matter


Particulate Matter = Extraneous particles, foreign particles

Risk of Harm = Probability of occurrence

× Detectability × Severity of harm



Aramus Fluoropolymer Bag vs. EVA Bag



- Entegris Aramus bag has the lowest cumulative particulate counts with particulates/mL $\geq 0.5 \mu\text{m}$
- Entegris Aramus SUS is 72X cleaner than EVA Product A, 38X cleaner than EVA Product B

CONCLUSION

- Entegris Life Sciences can leverage semiconductor purity and yield to biopharma/pharma manufacturing
- Entegris SUS products are manufactured in the cleanest environment
- Entegris' liquid particle counter AccuSizer SPOS system offers particulate monitoring down to 0.5 μm
- Aramus single-use systems has the lowest particulate counts comparing to market offerings of PE and EVA bags, especially down to 0.5 μm particle size



HOW TO MINIMIZE PARTICULATE CONTAMINATION

Quality by Design

- High-quality raw material with better particulate control
- Operator training for quality work
- In-process and lot release quality control (QC)
- ISO 9001 and ISO 13485 certification for quality management systems

Clean Manufacturing Environment



ISO 5 and ISO 7 clean room

Quality Control Testing



AccuSizer SPOS Systems, Entegris



Nicom Nano DLS/ZLS Systems, Entegris

FUTURE WORK

- In-line monitoring of particulate size for LNP production
- Profiling of particulate contamination in SUS assembly together with different tubings and connectors
- Subvisible particulate contamination's impact on purity, yield and potency of loaded viral vector
- Subvisible particulate contamination's impact on purity, yield and potency of loaded lipid nanoparticle